

THE MOST FLEXIBLE AND ADVANCED PRODUCTS IN IONTOPHORESIS

Trivarion's® unique, extremely flexible design conforms well to uneven treatment sites, and delivers ions with a precise pH balance up to 80 mA*min.




ADVANTAGES

- Thinnest, most flexible electrode on the market
- Maximum strength, precise buffering with new Carbon-Ag/AgCl technology
- Large ground electrode surface area lowers sensation and increases comfort
- Highly absorbent drug matrix – no messy hydrogel
- Large surface adhesion area




Item 001-48


TRIVARION® IS AVAILABLE IN THE FOLLOWING SIZES:



SMALL 1.5cc
Item 00115



MEDIUM 2.0cc
Item 00120



LARGE 3.5cc
Item 00135

**TRIVARION®
BUTTERFLY**



MEDIUM 2.0cc
Item 001-60

ACTIVADOSE® II CONTROLLER

- Two controls operate dose and current. Push-knob ON/OFF & PAUSE features. Improved display indicates dose, time and current. Dose range from 0 to 80 mA*min with 0.1 mA*min increments.
- Includes improved visual indicators and audible alerts.
- Time calculations are performed automatically and are recalculated if dose and current settings are changed.
- More power than dual-channel devices for better penetration; up to 80 volts.



Trivarion® and ActivaDose®II are registered trademarks of ActivaTek, Inc.

All of our Products are
**LATEX
FREE**



To learn more, please visit our website at
www.activatekinc.com

Phone & Fax: **1-800-680-5520**